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INSTRUCTION FOR USE

DE 2	EN 3	AR 4	BG 5	CS 6	DA 7	EL 8
ES 9	ET 10	FI 11	FR 12	HR 13	HU 14	IT 15
JA 16	LT 17	LV 18	NL 19	NO 20	PL 21	PT 22
RO 23	RU 24	SK 25	SL 26	SV 27	TR 28	ZH 29

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INTENDED USE

Cable Integrated Mounting systems (CIM) are designed to hold and position monitors, mouse/keyboard, laptops and other medical devices or appliances (hereinafter referred to as "ME devices"). Cables can be integrated inside the arms. The systems can be attached to wall, mobile carts, tables, and other surfaces or devices by specific mounting adapters. CIM systems support medical applications, but do not serve a medical purpose on its own. This instruction for use is intended for the operator and user of a CIM system. Instructions for ME devices and their accessories are not shown here.

Do not lean on CIM systems or ME equipment, as this may damage the equipment or cause instability. Risk of injury. Retain this instruction and observe all safety and other instructions.

Avoid penetration of disinfection and cleaning liquids into the interior of the system.

Do not use steel wool or abrasives.

The compatibility of disinfectant or cleaning fluids should be tested in advance on a non-critical spot.

Tests with the following groups of active substances have been carried out: alcohol-based rapid disinfectants (e.g. Bacillol AF/ Hartmann), aldehydes and quaternary ammonium compounds (e.g. B5/Orochemie), alcohol-free quaternary ammonium compounds (e.g. Sani-Cloth AF3/PDI Healthcare). The substances mentioned are VAH/ DGHM listed / EPA-registered and must be used according to the manufacturers' instructions.

The use of acetone and trichloroethylene will cause permanent damage of surfaces.

Cables integrated into the CIM system may also be affected by unsuitable agents. With any doubts, contact the manufacturer.

Never pull on cables or electrical lines. Damaged cables, lines or defect ME devices must be disconnected from the power supply immediately.

The swivel range is limited by internal rotational stops. Do not attempt to force beyond the swivel stop as this may cause damage to cables or equipment and endanger operators.

The swivel range of the CIM systems must be outside the patient area. No other devices may be located within swivel range of the arm. Risk of damage!

Move carefully and mind pinch points.

Stow all support arms in a safe position during transport.

Always move carts with both hands to avoid tipping.

Installed ME devices or cables can operate with dangerous electrical voltages.

Do not connect live cables to power until installation has been completed. Disconnect all live cables from power and equipment before carrying out inspection/disassembly.

Retain this instruction and observe all safety and other instructions.

CIM systems are maintenance-free.

Inspection should be performed by qualified personnel at least once a year.

Where necessary, CIM systems are equipped with locking levers/ parking brakes. To change the position, secure the ME device with your hand, release the locking lever, re-adjust the position and tighten the locking lever/ parking brake again.

Carts are equipped with one or more lockable castors. After positioning, lock all brakes by engaging the locking lever with your foot.

CIM systems are intended for use or storage in dry indoor locations.

Assembly, disassembly, cabling and installation of ME devices as well as inspection must be carried out by qualified personnel only. Corresponding assembly instructions must be followed.

It must be ensured that the chosen fixing or mounting adapter is suitable for the respective weight load. It is the customers/ system integrators responsibility to use adequate cables and ensure grounding/ earthing of the system.

CIM med GmbH recommends the consultation of your safety officer and biomedical technician to ensure safe installation and use of the CIM system.

CIM installation instructions including a inspection plan as well as instructions for use and be found here:

<https://cim-med.com/downloads>

Please also contact our service department for any questions.

**CE MARKING**

CIM systems and their accessories comply with applicable safety requirements.

**OBLIGATION TO REPORT**

All serious incidents occurring in connection with or caused by CIM systems must be reported to the manufacturer and the responsible authorities.

**MANUFACTURER**

CIM med GmbH

Margot-Kalinke-Str. 9, 80939 Munich, Germany

www.cim-med.com

**RECYCLING INFORMATION**

CIM systems contain valuable raw materials and can be recycled to a large extent. Electrical equipment must be collected and disposed of separately.

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